DESS

> DENTAL SMART SOLUTIONS

DESS

IFU

PROSTHETIC FIXTURES FOR DENTAL IMPLANTS

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DI-MUL-16-DSS

01-08-2019

Prosthetic Fixtures for dental implants

INSTRUCTIONS FOR USE

ENGLISH (INTERNATIONAL)

CAUTION – PLEASE READ CAREFULLY

DESS® products should only be used by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnosis, planning, dental surgery or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. DESS® products designed for single-use only must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health like patient tissue infection.

All DESS® components must be dry fitted before use to check the correct fitting. The clinician will be responsible for correct application of those restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor.

U.S. Federal law restricts those devices to sale by or on the order of a licensed dentist / physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must a screw check. If the screws have suffered unusual wear, the complete integrity of the implant-abutment should be checked. New screws must be used for any revision, correction or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

During any intraoral use and manipulation all DESS® products must be secured to prevent aspiration due to its shape and size.

INDICATIONS

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the dental laboratory.

CONTRAINDICATIONS

All materials used are biocompatible; however, some patients may present allergies or hypersensitivity to any of the materials and its components (specified in the table).

COMPATIBILITY INFORMATION

All DESS® components are available for a variety of connections. For compatibility with dental implants and implant analogues please follow our catalogue or guides or contact your local distributor.

To find necessary information about ANKYLOS® C/X system, please read the last paragraph of these instructions.

STERILIZATION AND RE-USE

All products are supplied NON STERILE. For sterilisation, we recommend autoclaving the product at 121°C for 30 minutes, drying time 30 minutes (in accordance with standard UNE-EN ISO 17665-1:2007). Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

ARTICLE/CLASS	GROUP DESS REF. N°	MATERIAL	STERILIZATION	RE-USE	
ANALOGUE / Class I	14.XXX and 14.XXX-P10	Surgical stainless steel AISI-303	In the autoclave	Not recommended	⊗
TRANSFER / Class I	17.XXX	Body: surgical stainless steel AISI-303 screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	After sterilization	-
SCREW / Class IIb	19.XXX and 19.XXX-P10	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CASTABLE ABUTMENT / Class IIb (28.XXX, 29.XXX), Class I (11.XXX-P10, 12.XXX-P10)	28.XXX, 29.XXX and 11.XXX-P10, 12.XXX-P10	Polyoxymethylene copolymer (POM)	N/A	Single use	⊗
GOLD UCLA / Class IIb	65.XXX and 66.XXX	Gold-Palladium alloy [Cendres+Métaux Ceramicor®]	Autoclave before use on the patient	Single use	⊗
Co-Cr BASE / Class IIb	30.XXX and 31.XXX	Cobalt-Chromium-Molybdenum alloy	Autoclave before use on the patient	Single use	⊗
INTERFACE* / Class IIb	15.XXX, 16.XXX, 36.XXX, 37.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CLINICAL ABUTMENTS* (STRAIGHT AND ANGLED)/ Class IIb	13.XXX, 22.XXX and 23.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
TRANSEPIOTHEIAL ABUTMENTS* (Multi-unit* and UniAbutment*) / Class IIb	40.XXX/X, 41.XXX/X, 43.XXX/X and 43.XXX/XR	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
HEALING ABUTMENT* / Class IIb	20.XXX/X and 20.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CONVERTER ABUTMENT / Class IIb	34.001 and 35.001	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
TEMPORARY ABUTMENT*/Class IIb	24.XXX and 25.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
SCAN ABUTMENT / Class I	50.XXX and 52.XXX	Polyetheretherketone (PEEK) For ANK System: Ti alloy ELI Ti-6Al-4V screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	In case of intraoral scanning: after sterilization	-
DESS LOC* / Class IIb	88.0XX/X and 88.XXX	Titanium alloy ELI Ti-6Al-4V, Zn coated	Autoclave before use on the patient	Single use	⊗
DESS LOC Retention Kit / Class I	DLRK-P1 and DLRK-P2, DLER-P1 and DLER-P2, DLN1-P4, DLN2-P4, DLN3-P4, DLN4-P4, DLN5-P4, DLN6-P4, DLN7-P4	Male – Polynil (Nylon) Cup – Titanium Alloy ELI Ti-6Al-4V Spacer ring – White Silicone	N/A	Single use	⊗
MASSIVE ABUTMENT / Class IIb	61.XXX and 71.XXX	Titanium alloy ELI Ti-6Al-4V (61.XXX) and Cobalt-Chromium-Molybdenum alloy (71.XXX)	Autoclave before use on the patient	Single use	⊗
TOOLS / Class I	MDTXXX and DTXXXXXX	Handles: surgical stainless steel SUS 316L Tips: stainless steel AISI 420 MOD	In the autoclave	In case of intraoral use: after sterilization	-

*Some parts are covered with MetAlive® coating and marked with ML letters on the beginning of the reference number.

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGUE

The analogues are for simulation of dental implant connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analogue are securely affixed. Check that the analogue and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analogue, because such a connection may be impossible to achieve with a product that may no longer be within specifications.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER

Transfer is used for closed and open tray technique, for the transfer of implant position from the intraoral situation to the model in dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove abutment, clean the connection with spray water/air and dry. Prepare transfer adequate to implant connection and platform size. Place it on implant and check correct position. Put the screw and fix manually.

- For closed tray technique choose the short screw. After tightening block the hex of the screw with wax. Use hydrocolloid, polyethylene or soft silicone. After hardening place the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For open tray technique choose the long screw. After tightening place the impression tray to ensure full access to transfer screws from outside. Before impression taking the transfers can be ferulized with acrylic resin placed over dental floss between consecutive implants. After hardening release all screws and take out the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over implant or analogue. Make sure to secure the parts with corresponding screw and observe specified torque value placed on the label or ask your local DESS® distributor. You can download torques brochure from <http://www.dess-abutments.com/downloads>. For best results, the following conditions must be meticulously met:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The driver should be placed in the longitudinal axle of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every checks thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case.
- **Position the patient to avoid aspiration in case the screw falls during screwing/unscrewing.**
- Check compatibility of the screw with the implant model to which it will be connected.

Some screw models are available with DLC coating [Diamond-Like Carbon, Chromium Carbonitride coating – CrCN] to provide a low friction surface, which improves the preload of the screw and provide better hold. The Chromium Carbonitride coating is not toxic and has excellent mechanical, tribological and biological properties. DLC is often used in blood contacting implants such as stents and heart valves and to reduce wear in load bearing joints.

APPLICATIONS AND DIRECTIONS FOR USE – CASTABLE ABUTMENT

The castable abutment is for prosthetic restorations prepared by dental technicians in the dental laboratory. For greater connection precision use metal abutments which have been previously mechanised. However, if you decide to use the casting method, please take the following precautions:

- Use a flash of wax that will create sufficient gap around the abutment to compensate for the greater dilatation coefficient between the castable material and its surroundings.
- Gentle torque to prevent deformation.
- Perform castings of a size and form that favours filling of the cavities and prevents the appearance of air bubbles.
- Use materials with high fluidity for castings in conflictive models

APPLICATIONS AND DIRECTIONS FOR USE – GOLD UCLA

The Gold UCLA is used for prosthetic restorations prepared by dental technicians in the dental laboratory. The gold-plastic abutment is an implant abutment that consists of a cast-on gold alloy base and a fully burn-out plastic sleeve. The laboratory processing of the gold-plastic abutment utilizes the cast-on technique. During casting, the gold portion is suffused on the proposed contact areas by melting of the cast-on alloy. Metallic joint is obtained through diffusion. Depending of the model differ weight of pre-milled part (see table). **Processing:** The abutment must be carefully hand-tightened into the lab analogue with a lab screw. The lab screw is designed exclusively for laboratory processing of the abutment. The existing abutment screw should be used only for the final integration of the restoration after being sterilized. Before placing casting sprues the metal base should be cleaned with a cotton swab and alcohol. Casting sprues should be placed possibly in the long axis of the crown to avoid the bubble concentration inside the construction. The elements need to be placed away from casting ring heating center. While setting the sprue, ensure correct positioning of the wax-up casting object in the casting ring. The metal tube of the abutment should not be placed too close to the wall or floor of the ring to prevent heat removal (cooling vents). If the metal tube cools too rapidly, this can lead to defective casts. The use of speed investment is not recommended. The investment expansion should be adjusted to zero expansion. Make sure that the investment is free of air bubbles.

Observe the instructions of the investment material's manufacturer. Since the abutment can be exposed to very high loading in the patient's mouth, only high-strength, precious metal alloys should be used (casting, universal, or metal ceramic alloys). To ensure full casting, the final temperature during preheating of the casting muffle must be at least 750 °C for casting and universal alloys and at least 900 °C for metal ceramic alloys. The final temperature should be held for 30 to 60 minutes according to the size of the casting ring. The casting temperature of high precious metal alloys must be at least 150 °C below the solidus temperature of the gold-plastic abutment (1400 °C). This prevents melting of the gold-plastic abutment.

Observe the instructions of the alloy manufacturer. To prevent stresses in the cast structure resulting from too rapid cooling, cool the casting muffle at room temperature. It is better to deflask the casting mould using an acid bath in an ultra sound cleaner.

Sandblasting can damage the settlement. During neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. If working with cemented restorations, you should develop the crown and bridge frameworks on the abutments in the same way as for periodontic restorations.

COMPOSITION (percentage by mass)

Au 60.0 %
Pd 20.00 %
Pt 19.00 %
Ir 1.00 %

TECHNICAL DATA

Yield strength Rp 0.2 > 650 MPa
Tensile strength Rm > 750 MPa
Hardness HV5 > 215
Fracture elongation > 2 %
Melting range 1400 - 1490 °C
CTE (25 – 500 °C) 11.9 x 10⁻⁶ K⁻¹
Color white

ARTICLE

DESS REF. N°

BASE WEIGHT (g)

Octagon Gold UCLA Non-Engaging on Implant	65.009	0.303
Octagon Gold UCLA Engaging on Implant	66.009	0.374
Tri-Lobe Gold UCLA Engaging NP	66.004	0.246
Tri-Lobe Gold UCLA Engaging RP	66.005	0.326
Internal Hex USA Gold UCLA Engaging NP	66.017	0.372
Internal Hex USA Gold UCLA Engaging RP	66.018	0.334
Tri-Lobe Gold UCLA Non-Engaging NP	65.004	0.152
Tri-Lobe Gold UCLA Non-Engaging RP	65.005	0.182
Internal Hex USA Gold UCLA Non-Engaging NP	65.017	0.325
Internal Hex USA Gold UCLA Non-Engaging RP	65.018	0.287

APPLICATIONS AND DIRECTIONS FOR USE – Co-Cr BASE

The Co-Cr Base is used for prosthetic restorations prepared by dental technicians in the dental laboratory. The Co-Cr abutment is an implant abutment that consists of a cast-on Co-Cr alloy base and a fully burn-out plastic sleeve. The laboratory processing of the Co-Cr base abutment utilizes the cast-on technique. During casting, the Co-Cr portion is suffused on the proposed contact areas by melting of the cast-on alloy. Metallic joint is obtained thanks to the retention of the grooves of the Co-Cr base and the suffused metal. **Processing:** The abutment must be carefully hand-tightened into the lab analogue with a lab screw. The lab screw is designed exclusively for laboratory processing of the abutment. The existing abutment screw should be used only for the final integration of the restoration after being sterilized. Create a wax-collar of at least 0.3 mm thickness above the metal indexing feature of the base. Before placing casting sprues the metal base should be cleaned with a cotton swab and alcohol. Casting sprues should be placed possibly in the long axis of the crown to avoid the bubble concentration inside the construction. The investment should freely pass through the screw channel. The elements need to be placed away from casting ring heating center. While setting the sprue, ensure correct positioning of the wax-up casting object in the casting ring. The metal tube of the abutment should not be placed too close to the wall or floor of the ring to prevent heat removal (cooling vents). The use of speed investment is not recommended. The investment expansion should be adjusted to zero expansion. Make sure that the investment is free of air bubbles. Heat parameters need to be tuned as for conventional CoCr casting to avoid any imperfections in the construction. Too rapid cooling of the metal tube can lead to defective casts. Observe the instructions of the investment material manufacturer. To ensure full casting, the final temperature during preheating of the casting muffle must be at least 950 °C. The final temperature should be held for 30 to 60 minutes according to the size of the casting ring. To prevent stresses in the cast structure resulting from too rapid cooling, cool the casting muffle at room temperature. It is better to deflask the casting mould using an acid bath in an ultra sound cleaner. Sandblasting can damage the settlement. During neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. If working with cemented restorations, you should develop crown and bridge frameworks on the abutments in the same way as for periodontic restorations.

COMPOSITION (percentage by mass)

Co 65.4 %
Cr 27.75 %
Mo 5.06 %
Additional elements less than 1 % (Fe, Mn, Si, N, Ni, C, Ti, P, S)

TECHNICAL DATA

Yield strength¹ Rp 0.2 > 827 MPa
Tensile strength¹ Rm > 1172 MPa
Hardness¹ HV10 459
Fracture elongation¹ A5 21 %
Solidus¹ 1369 °C
Liquidus¹ 1415 °C
CTE¹ (25 – 500 °C) 14.1 x 10⁻⁶ K⁻¹
Color white

¹ Variations are possible depending on the batch.

APPLICATIONS AND DIRECTIONS FOR USE – INTERFACE

The interface is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium interface is to support the bridge or individual zirconium dioxide elements made in CAD/CAM or manual milling machine. It can also be used as a connector between full-shape screwed crowns (from any material) and an implant, either as narrow abutment in some cases. Every titanium interface comes with a calibrated POM castable which can be burned out without residues if necessary. **Contraindications:** do not use the titanium interface for restorations with cantilever on a single implant, bruxism, insufficient space or direct metal-to-interface casting.

To use this product in digital dentistry DESS® library and scan abutment are needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Some dental CAD/CAM systems require a double scan with diagnostic wax-up. If you need additional assistance with installation and/or proper use of the files and with use of piece itself, please contact your DESS® local distributor or find necessary information in corresponding instructions. All instructions can be downloaded from our website <http://www.dess-abutments.com/downloads>. To improve cement adhesion we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following recommendations of cement manufacturer. The ceramic surface of mesostructure in cementing zone should be sandblasted and cleaned/degreased. For secure grip the diameter of interface and its height should not be reduced e.g. by grinding.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL ABUTMENTS STRAIGHT AND ANGLED

The straight and angled abutments are used for the prosthetic restorations prepared by dental technicians in dental laboratory or by dentists at chairside in dental clinic. Straight and angled abutments form the core of prosthetic crowns and bridges or are used for their support. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts and that they are in the correct place and position for all secondary parts. We strongly recommend a radiography to verify the correct settlement after tightening. **Caution:** If customization is needed, do not grind the wall to less than 0.5 mm thick. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before cementing the final restoration, check the torque according to implant connection and size. When cementing use definitive luting cement or use temporary cement after sealing of access hole with Teflon tape or gutta-percha. Remove cement excess from edge of the crown to avoid peri-implantitis which can lead to implant loss. After cementation remember to check static and dynamic occlusion using conventional protocol.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSEPITHELIAL ABUTMENTS: Multi-unit® and UniAbutment®

Those abutments are used for prosthetic restorations prepared by dentist in dental clinic. A pre-manufactured dental implant abutment directly connected to the dental implant intended for use as an aid in prosthetic rehabilitation for multiple fixed-removable prosthesis or screw-retained. We strongly recommend to verify the final abutment seating using radiographic imaging.

DESS® offers two systems for this kind of rehabilitation: Multi-unit® system and UniAbutment®. In case of straight Multi-unit® abutment: place appropriate abutment. Use plastic holder to facilitate the insertion. Tighten the abutment using DESS® key for Multi-unit® and manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). In case of angulated Multi-unit® abutment: place appropriate abutment using metal holder at desired position and tighten the abutment using system specific screwdriver in manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). Take an impression in a standard procedure and provisionalize the patient. If temporary prosthesis is not necessary, place healing caps. In case of UniAbutment®: place adequate abutment choosing different type (20° or 45° - depends of implant divergence) and height. Use attached holder screw to facilitate the insertion. Tighten the abutment using DESS® key for UniAbutment® and take impression using standard techniques (for recommended torque use Clinical Screw Recommended Torque document).

APPLICATIONS AND DIRECTIONS FOR USE – HEALING ABUTMENT

The healing abutment is used on the patient as an auxiliary transmucosal abutment. It is placed on the implant fixture before prosthetic restoration to facilitate the formation of a soft tissue sulcus. Choose the correct height to ensure its correct function and prevent the transmission of chewing forces. Before placing ensure that implant platform stays free from any tissue residues. Use gentle manual torque when installing.

APPLICATIONS AND DIRECTIONS FOR USE – CONVERTER ABUTMENT

The converter abutment is used in prosthetic restorations prepared by dentist in dental clinic or by technicians in dental laboratory. Converter abutments are used for changing the connection from internal to external through the Morse taper connection to secure implant-to-abutment connection. These elements are placed over the dental implant acting as part of the prosthesis and are connected to the implant by a captive screw.

After placing in implant, use adequate DESS® part in ratchet to tighten the abutment: for non-engaging abutments use octagonal key. For engaging abutments use corresponding screwdriver tip in adapter.

For both adapters torque recommended can be found in Clinical Screw Recommended Torque document, continue with standard prosthetic procedure.

APPLICATIONS AND DIRECTIONS FOR USE – TEMPORARY ABUTMENT

The temporary abutment is used for prosthetic restorations prepared by dentist at the chairside, works as support for provisional screwed crowns, bridges and prosthesis. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend a radiography to verify the correct settlement after tightening. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before placing the temporary restoration check the torque according to implant connection and size. It's important to seal the access hole with Teflon tape or gutta-percha and composite. Any material excess near implant platform should be removed, as it can provoke peri-implantitis and can lead to implant loss. It is very important to check static and dynamic occlusion avoiding excessive loading forces which may prolong or prevent the proper osseointegration, especially in cases of immediate load.

APPLICATIONS AND DIRECTIONS FOR USE – SCAN ABUTMENT

Scan abutment is used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentist at the chairside in the CAD/CAM process. Indicated for obtain geometric data from the master model using a desktop laboratory 3D scanner or for optical impressions using intraoral 3D scanner. For proper process a DESS® digital library is needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Before scanning the abutment and analogue should be visually checked for surface damage or implant connection for any tissue residues. For greater scanning precision, we recommend to locate the flat surface of the scan abutment in palatal/lingual orientation. Fasten the abutment using corresponding screw (indications in our catalogue and on the product label) by hand or with maximum 10 Ncm. In case of intraoral scan abutment, fix it with fingers using embedded screw. Scan abutment is a precision tool and being tightened hard may change its geometry causing errors in scanning process and discrepancy in accuracy. For most of scanners a scan-enhancing spray is not required.

The scanning process should be performed as recommended by CAD/CAM system manufacturer. It is important to choose the correct implant connection in software and corresponding type for the chosen restoration (engaging/non-engaging). After process the scan abutment can be loosened and placed gently on the tray or box. To complete the scanning gathering additional information (e.g. occlusion index, silicone bite, gingiva shape) can be necessary. In case of intraoral use a sterilization is needed.

APPLICATIONS AND DIRECTIONS FOR USE – DESS LOC

DESS LOC is a resilient implant attachment used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentist at the chairside. Indicated for removable dentures.

Procedure: Measure the tissue thickness from the implant platform to the highest side of the gingival crest near the implant. Choose the corresponding height or 1.5 mm higher (the attachment should not be submerged) and place using recommended torque on Clinical Screw Recommended Torque document. The retention male part can be mounted in dental clinic as well as in dental laboratory with self-curing resin using standard techniques.

APPLICATIONS AND DIRECTIONS FOR USE – DESS LOC RETENTION KIT

DESS LOC retention system is a part of the anchoring system for full or partial overdentures on implants, in cases where conventional full dentures have problems with retention and stability. The male retention kit can be mounted in the dental clinic using the self-curing resin, and in the dental laboratory, using the standard procedure. **Direct method:** Place the silicone spacer over the attachment covering the entire area around. If there is any space between the cup and soft tissue, put more spacers until it disappears. Place the cup with the black Nylon male on the DESS LOC Attachment in patient's mouth. Prepare a site for the cup in the prosthesis and perforate the buccal side to evacuate the excess of material. Apply the dual- or cold-curing resin around the cup and within the holes. Place the prosthesis in place while keeping the height of the occlusion and hold in position until the resin hardens. Remove the denture, spacers and cut of the excess. Polish before changing to the final male Nylon. **Indirect method:** same as direct, but on the master model with 14.054 analogs. **Retention force/divergence between implants:** blue (1.5 lbs/680g, up to 20°), rose (3 lbs/1360g, up to 20°), transparent (5 lbs/2270g, up to 20°), black (Processing Nylon male, 0 lbs, 0g, up to 40°), red (1 lbs/450g, up to 40°), orange (2 lbs/910g, up to 40°), green (4 lbs/1810g, up to 40°). DLRK-P1 and DLRK-P2 are kits up to 20°, DLER-P1 and DLER-P2 are kits up to 40°.

APPLICATIONS AND DIRECTIONS FOR USE – METALIVE® COATING

Some DESS® parts are covered with MetAlive® coating which enables protein-mediated attachment of the gingival tissue to the abutment surface. However, this does not change the general indications for use of the modified abutment.

Warning: Do not treat the gingival contact surface with any bonding agents, cements, cleaning agents or disinfectants. Do not clean the surface using abrasive technologies. Plaque removal should be done using soft scaling instruments. If sandblasted in laboratory, the gingival contact area must be protected with wax. Before sterilization or for cleaning, water vapour can be used as well as acetone and absolute alcohol in ultrasonic cleaner for 5 min. After cleaning or sterilizing do not touch surface with bare fingers and keep the surface hygienic. Please avoid disinfectants such as iodine, Hydrogen peroxide, chlorhexidine, polyvidone, etc. Avoid pointing strong laser or IR-radiation to modified surface. Do not remove MetAlive® surface material from gingival contact area with dental burr or by polishing instruments such as diamond rubber. Be sure that the surface is not in contact with hard sharp edges or points. **Processing:** MetAlive® abutments can be sterilized in autoclave at 121° C. After sterilization have the gingival tissue wound refreshed to enhance tissue repair and attachment, also after replacement of healing abutment with the permanent MetAlive® abutment. Ensure maximal initial contact with gingival tissue. Ensure minimal movement of gingival tissue during healing. Patient should use soft toothbrushes and non-abrasive toothpaste. For cleaning the interdental spaces dental floss or interdental brushes can be used.

APPLICATIONS AND DIRECTIONS FOR USE – MASSIVE ABUTMENT

Massive abutment is used for prosthetic restorations prepared by dental technicians in dental laboratory. The pre-milled implant-abutment interface of the blank is fabricated on an industrial lathe in order to provide consistent and accurate results. The connection is more accurate and the abutment shape can be milled with burs to meet specific patient requirements.

Caution: For proper milling a DESS fixture adapter is needed. The minimum wall thickness shall not be less than 0.4 mm.

APPLICATIONS AND DIRECTIONS FOR USE – TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. By use of exchangeable screwdriver tips may be used for various implant systems. For intraoral use, all tools must be thoroughly autoclaved before use. Where using removable heads, make sure that it cannot come loose from the driver. The driver should always be fastened with a safety pin and cord attached the groove designed for this purpose. It should be handled with care, taking the necessary precautions not to hurt the patient. Before using the driver, check that it is the correct size and shape for the screw head that will be used. Only use drivers that are in perfect conditions. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

In case of adaptors please assure that corresponds to 4 x 4 connection and to its abutment system (Multi-unit®, UniAbutment®, octagonal or Locator®).

INTERNAL ANK SPECIAL REQUIREMENTS

The thickness of the anti-rotation elements is reduced due to the widened screw channel in order to accommodate the use of a normal screw. To prevent bending of those parts under applied pressure, please ensure that the screw is threaded through the abutment whilst restoration takes place and before placing it in the analog or implant. If this precaution is taken, the anti-rotation elements remain supported and safe from shear and compression forces. If, for any reason the screw needs to be removed, place it back again in its position before applying any loads on the abutment.

EXPLANATION OF THE SYMBOLS THAT APPEAR ON THE LABELS

	For single use only		Batch No.		Reference
	NON Sterile product		Read instructions	 0051	Medical product compliance marking controlled by the notified body number 0051 (class IIb)
Rx Only			Use only on prescription		
	Manufacturer information		Manufacturer: Terrats Medical SL Calle Mogoda 75-99 08210 – Barberá del Vallés -Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 info@dess-abutments.com		

All DESS® products are manufactured according to ISO 9001 and ISO 13485 and Class IIb bare the CE marking.

3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA
 ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE® and FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants
 BioHorizons® registered trademark of BioHorizons, INC.
 Brånemark System®, NobelReplace® and Replace Select™ and NobelActive™ and NobelReplace® CC and NobelSpeedy™ and Multi-unit® registered trademarks of Nobel Biocare AB, Sweden
 ITI®, Straumann® and synOcta® registered trademarks of Straumann Holding AG, Switzerland
 MIS® registered trademark of MIS Implant Technologies Ltd
 Zimmer Screw-Vent® and SwissPlus® registered trademarks of Zimmer Dental INC, USA
 Dyna® registered trademark of Dyna Dental Engineering BV
 Locator® registered trademark of Zest Anchors INC.

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. Some DESS® items may not be available in your country. These instructions replace all previous editions. For any detailed information about DESS® products please contact your local distributor.

Prosthetic Fixtures for dental implants

INSTRUCTIONS FOR USE

ENGLISH (CANADA)

WARNING – PLEASE READ CAREFULLY

DESS® products are intended to be used only by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnosis, planning, dental surgery or prosthetic techniques. If in doubt regarding the device intended use, please contact the manufacturer. DESS® products, designed for single-use only, must never be re-used. If reused, there is a risk of damage to the device or deterioration of its characteristics that could lead to failure of the prosthetic solution failure and / or to deterioration of the patient's health such as patient tissue infection.

All DESS® components must be dry fitted before use to check the correct fitting. The clinician is responsible for the correct application of dental restorative devices as well as for planning and conducting those dental procedures that are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please call the distributor.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must include a screw check. If the screws are observed to have suffered unusual wear, the complete integrity of the implant-abutment should be checked. New screws must be used for any revision, correction or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

During any intraoral use and manipulation, all DESS® products must be secured to prevent aspiration due to its shape and size.

INDICATIONS

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the laboratory.

CONTRAINDICATIONS

All materials used are biocompatible; however some patients may present allergies or hypersensitivity to any of the materials and its components (specified in the table).

COMPATIBILITY INFORMATION

All DESS® components are available for a variety of connections. For compatibility with dental implants and implant analogues please follow our catalogue or contact your local distributor. To find necessary information about ANKYLOS® C/X system, please read the last paragraph of these instructions.

STERILIZATION AND RE-USE

All products are supplied NON STERILE. For sterilisation, we recommend autoclaving the product at 121°C for 30 minutes, drying time 30 minutes (in accordance with standard UNE-EN ISO 17665-1:2007). Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

ARTICLE/CLASS#	GROUP DESS REF. N°	MATERIAL	STERILIZATION	RE-USE	
ANALOGUE / Class I	14.XXX and 14.XXX-P10	Surgical stainless steel AISI-303	In the autoclave	Not recommended	⊗
TRANSFER / Class I	17.XXX	Body: surgical stainless steel AISI-303 screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	After sterilization	-
SCREW / Class III	19.XXX and 19.XXX-P10	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CASTABLE ABUTMENT / Class III (28.XXX, 29.XXX), Class I (11.XXX-P10, 12.XXX-P10)	28.XXX, 29.XXX and 11.XXX-P10, 12.XXX-P10	Polyoxymethylene copolymer (POM)	N/A	Single use	⊗
GOLD UCLA / Class III	65.XXX and 66.XXX	Gold-Palladium alloy [Cendres+Métaux Ceramicor®]	Autoclave before use on the patient	Single use	⊗
Co-Cr BASE / Class III	30.XXX and 31.XXX	Cobalt-Chromium-Molybdenum alloy	Autoclave before use on the patient	Single use	⊗
INTERFACE* / Class III	15.XXX, 16.XXX, 36.XXX, 37.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CLINICAL ABUTMENTS* (STRAIGHT AND ANGLED) / Class III	13.XXX, 22.XXX and 23.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
TRANSEPIHELIAL ABUTMENTS* (Multi-unit® and UniAbutment®) / Class III	40.XXX/X, 41.XXX/X, 43.XXX/X and 43.XXX/XR	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
HEALING ABUTMENT* / Class III	20.XXX/X and 20.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
CONVERTER ABUTMENT / Class III	34.001 and 35.001	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
TEMPORARY ABUTMENT* / Class III	24.XXX and 25.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
SCAN ABUTMENT / Class I	50.XXX and 52.XXX	Polyetheretherketone (PEEK) For ANK System: Ti alloy ELI Ti-6Al-4V screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	In case of intraoral scanning: after sterilization	-
DESS LOC* / Class III	88.0XX/X and 88.XXX	Titanium alloy ELI Ti-6Al-4V, ZnN coated	Autoclave before use on the patient	Single use	⊗
PRE-MILLED ABUTMENT / Class III	61.XXX and 71.XXX	Titanium alloy ELI Ti-6Al-4V (61.XXX) and Cobalt-Chromium-Molybdenum alloy (71.XXX)	Autoclave before use on the patient	Single use	⊗
TOOLS / Class I	MDTXXX and DTXXXXXX	Handles: surgical stainless steel SUS 316L Screwdriver: stainless steel AISI 420 MOD	In the autoclave	In case of intraoral use: after sterilization	-

*Classification in accordance with Canadian Medical Devices Regulations

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGUE

The analogues are for simulation of dental implant connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analogue are securely affixed. Check that the analogue and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analogue, because such a connection may be impossible to achieve with a product that may no longer be within specifications.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER

Transfer is used for closed and open tray technique, for the transfer of implant position from the intraoral situation to the model in dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove abutment, clean the connection with spray water/air and dry. Prepare transfer adequate to implant connection and platform size. Place it on implant and check correct position. Put the screw and fix manually.

- For closed tray technique choose the short screw. After tightening, block the hex of the screw with wax. Use hydrocolloid, polyethylene or soft silicone. After hardening, place the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For open tray technique choose the long screw. After tightening, place the impression tray to ensure full access to transfer screws from outside. Before impression taking, the transfers can be ferulized with acrylic resin placed over dental floss between consecutive implants. After hardening release all screws and take out the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over an implant or analogue. Make sure to secure the parts with corresponding screw and observe the specified torque value indicated on the label or ask your local DESS® distributor. You can download torques brochure from <http://www.dess-abutments.com/downloads>. For best results, the following conditions must be meticulously met:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The driver should be placed in the longitudinal axle of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every check thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case. If you have any doubts about the fit of the screw, consult the local distributor.
- Position the patient to avoid aspiration in case the screw falls during screwing/unscrewing.**
- Check compatibility of the screw with the implant model to which it will be connected.

Some screw models are available with DLC coating [Diamond-Like Carbon, Chromium Carbonitride coating – CrCN] to provide a low friction surface, which improves the preload of the screw and provide better hold. The Chromium Carbonitride coating is not toxic and has excellent mechanical, tribological and biological properties. DLC is often used in blood contacting implants such as stents and heart valves and to reduce wear in load bearing joints.

APPLICATIONS AND DIRECTIONS FOR USE – CASTABLE ABUTMENT

The castable abutment is for prosthetic restorations prepared by dental technicians in the dental laboratory. For greater connection precision use metal abutments which have been previously mechanised. However, if you decide to use the casting method, please take the following precautions:

- Use a flash of wax that will create sufficient gap around the abutment to compensate for the greater dilatation coefficient between the castable material and its surroundings.
- Use gentle torque to prevent deformation.
- Perform castings of a size and form that favour filling of the cavities and prevents the appearance of air bubbles.
- Use materials with high fluidity for castings in conflictive models

APPLICATIONS AND DIRECTIONS FOR USE – GOLD UCLA

The Gold UCLA is used for prosthetic restorations prepared by dental technicians in the dental laboratory. The gold-plastic abutment is an implant abutment that consists of a cast-on gold alloy and a fully burn-out plastic sleeve. The laboratory processing of the gold-plastic abutment utilizes the cast-on technique. During casting, the gold portion is suffused on the proposed contact areas by melting of the cast-on alloy. Metallic joint is obtained through diffusion. Depending on the model the weight of the pre-milled part differs (see table). **Processing:** The abutment must be carefully hand-tightened into the lab analogue with a lab screw. The lab screw is designed exclusively for laboratory processing of the abutment. The existing abutment screw should be used only for the final integration of the restoration after being sterilized. Before placing casting sprues the metal base should be cleaned with a cotton swab and alcohol. Casting sprues should be placed possibly in the long axis of the crown to avoid the bubble concentration inside the construction. The investment should freely pass through the screw channel. The elements need to be placed away from casting ring heating center. While setting the sprue, ensure correct positioning of the wax-up casting object in the casting ring. The metal tube of the abutment should not be placed too close to the wall or floor of the ring to prevent heat removal (cooling vents). If the metal tube cools too rapidly, this can lead to defective casts. The use of speed investment is not recommended. The investment expansion should be adjusted to zero expansion. Make sure that the investment is free of air bubbles. Observe the instructions of the investment material's manufacturer. Since the abutment can be exposed to very high loading in the patient's mouth, only high-strength, precious metal alloys should be used (casting, universal, or metal ceramic alloys). To ensure full casting, the final temperature during preheating of

COMPOSITION (percentage by mass)	TECHNICAL DATA
Au 60.0 % Pd 20.00 % Pt 19.00 % Ir 1.00 %	Yield strength Rp 0.2 > 650 MPa Tensile strength Rm > 750 MPa Hardness HV5 > 215 Fracture elongation > 2 % Melting range 1400 - 1490 °C CTE (25 – 500 °C) 11.9 x 10 ⁻⁶ K ⁻¹ Color white

ARTICLE	DESS REF. N°	BASE WEIGHT (g)
Octagon Gold UCLA Non-Engaging on Implant	65.009	0.303
Octagon Gold UCLA Engaging on Implant	66.009	0.374
Tri-Lobe Gold UCLA Engaging NP	66.004	0.246
Tri-Lobe Gold UCLA Engaging RP	66.005	0.326
Internal Hex USA Gold UCLA Engaging NP	66.017	0.372
Internal Hex USA Gold UCLA Engaging RP	66.018	0.334
Tri-Lobe Gold UCLA Non-Engaging NP	65.004	0.152
Tri-Lobe Gold UCLA Non-Engaging RP	65.005	0.182
Internal Hex USA Gold UCLA Non-Engaging NP	65.017	0.325
Internal Hex USA Gold UCLA Non-Engaging RP	65.018	0.287

the casting muffle must be at least 750 °C for casting and universal alloys and at least 900 °C for metal ceramic alloys. The final temperature should be held for 30 to 60 minutes according to the size of the casting ring. The casting temperature of high precious metal alloys must be at least 150 °C below the solidus temperature of the gold-plastic abutment (1400 °C). This prevents melting of the gold-plastic abutment. Observe the instructions of the alloy manufacturer. To prevent stresses in the cast structure resulting from cooling too quickly, cool the casting muffle at room temperature. It is better to deflask the casting mould using an acid bath with cemented restorations, you should develop the crown and bridge frameworks on the abutments in the same way as for periodontic restorations.

APPLICATIONS AND DIRECTIONS FOR USE – Co-Cr BASE

The Co-Cr Base is used for prosthetic restorations prepared by dental technicians in the dental laboratory. The Co-Cr abutment is an implant abutment that consists of a cast-on Co-Cr alloy base and a fully burn-out plastic sleeve. The laboratory processing of the Co-Cr base abutment utilizes the cast-on technique. During casting, the Co-Cr portion is suffused on the proposed contact areas by melting of the cast-on alloy. Metallic joint is obtained thanks to the retention of the grooves of the Co-Cr base and the suffused metal. **Processing:** The abutment must be carefully hand-tightened into the lab analogue with a lab screw. The lab screw is designed exclusively for laboratory processing of the abutment. The existing abutment screw should be used only for the final integration of the restoration after being sterilized. Create a wax-collar of at least 0.3 mm thickness above the metal indexing feature of the base. This prevents cracks in ceramic layer. Before placing casting sprues the metal base should be cleaned with a cotton swab and alcohol. Casting sprues should be placed possibly in the long axis of the crown to avoid the bubble concentration inside the construction. The investment should freely pass through the screw channel. The elements need to be placed away from casting ring heating center. While setting the sprue, ensure correct positioning of the wax-up casting object in the casting ring. The metal tube of the abutment should not be placed too close to the wall or floor of the ring to prevent heat removal (cooling vents). The use of speed investment is not recommended. The investment expansion should be adjusted to zero expansion. Make sure that the investment is free of air bubbles. Heat parameters need to be adjusted as for conventional CoCr casting to avoid any imperfections in the construction. Cooling the metal tube too quickly can lead to defective casts. Observe the instructions of the investment material manufacturer. To ensure full casting, the final temperature during preheating of the casting muffle must be at least 950 °C. The final temperature should be held for 30 to 60 minutes according to the size of the casting ring. To prevent stresses in the cast structure resulting from too rapid cooling, cool the casting muffle at room temperature. It is better to deflask the casting mould using an acid bath in an ultra sound cleaner. Sandblasting can damage the settlement. During neck polishing the connection surface needs to be covered with the analogue: this prevents any damage to the settlement part. If working with cemented restorations, you should develop crown and bridge frameworks on the abutments in the same way as for periodontic restorations.

COMPOSITION (percentage by mass)	TECHNICAL DATA
Co 65.4 % Cr 27.75 % Mo 5.06 % Additional elements less than 1 % (Fe, Mn, Si, N, Ni, C, Ti, P, S)	Yield strength ¹ Rp 0.2 > 827 MPa Tensile strength ¹ Rm > 1172 MPa Hardness ¹ HV10 459 Fracture elongation ¹ A5 21 % Solidus ¹ 1369 °C Liquidus ¹ 1415 °C CTE ¹ (25 – 500 °C) 14.1 x 10 ⁻⁶ K ⁻¹ Color white

¹ Variations are possible depending on the batch.

APPLICATIONS AND DIRECTIONS FOR USE – INTERFACE

The interface is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium interface is to support the bridge or individual zirconium dioxide elements made in CAD/CAM or manual milling machine. It can also be used as a connector between full-shape screwed crowns (from any material) and an implant, either as narrow abutment in some cases. Every titanium interface comes with a calibrated POM castable which can be burned out without residues if necessary. **Contraindications:** do not use the titanium interface for restorations with cantilever on a single implant, bruxism, insufficient space or direct metal-to-interface casting.

To use this product in digital dentistry, DESS® library and scan abutments are needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Some dental CAD/CAM systems require a double scan with diagnostical wax-up. If you need additional assistance with installation and/or proper use of the files and with use of piece itself, please contact your DESS® local distributor or find necessary information in corresponding instructions. All instructions can be downloaded from our website <http://www.dess-abutments.com/downloads> To improve cement adhesion we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following recommendations of cement manufacturer. The ceramic surface of mesostructure in cementing zone should be sandblasted and cleaned/degreased. For secure grip the diameter of interface and its height should not be reduced e.g. by grinding.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL ABUTMENTS STRAIGHT AND ANGLED

The straight and angled abutments are used for the prosthetic restorations prepared by dental technicians in dental laboratory or by dentists at chairside in dental clinic. Straight and angled abutments form the core of prosthetic crowns and bridges or are used for their support. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts and that they are in the correct place and position for all secondary parts. We strongly recommend an X-ray to verify the correct settlement after tightening. **Caution:** If customization is needed, do not grind the wall to less than 0.5 mm thick. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Position the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before cementing the final restoration, check the torque according to implant connection and size. When cementing use definitive luting cement or use temporary cement after sealing of access hole with Teflon tape or gutta-percha. Remove cement excess from edge of the crown to avoid peri-implantitis which can lead to implant loss. After cementation remember to check static and dynamic occlusion using conventional protocol.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSEPITHELIAL ABUTMENTS : Multi-unit® and UniAbutment®

Those abutments are used for prosthetic restorations prepared by dentist in dental clinic. A pre-manufactured dental implant abutment directly connected to the dental implant intended for use as an aid in prosthetic rehabilitation for multiple fixed-removable prosthesis or screw-retained. We strongly recommend to verify the final abutment seating using radiographic imaging.

DESS® offer two systems for this kind of rehabilitation: Multi-unit® system and UniAbutment®. In case of straight Multi-unit® abutment: place appropriate abutment. Use plastic holder to facilitate the insertion. Tighten the abutment using DESS® key for Multi-unit® and manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). In case of angulated Multi-unit® abutment: place appropriate abutment using metal holder at desired position and tighten the abutment -using system specific screwdriver in manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). Take an impression in a standard procedure and provisionalize the patient. If temporary prosthesis is not necessary, place healing caps. In case of UniAbutment®: place adequate abutment choosing different type (20° or 45° - depends of implant divergence) and height. Use attached holder screw to facilitate the insertion. Tighten the abutment using DESS® key for UniAbutment® and take impression using standard techniques (for recommended torque use Clinical Screw Recommended Torque document).

APPLICATIONS AND DIRECTIONS FOR USE – HEALING ABUTMENT

The healing abutment is used on the patient as an auxiliary transmucosal abutment. It is placed on the implant fixture before prosthetic restoration to facilitate the formation of a soft tissue sulcus. Choose the correct height to ensure its correct function and prevent the transmission of chewing forces. Before placing ensure that implant platform stays free from any tissue residues. Use gentle manual torque when installing.

APPLICATIONS AND DIRECTIONS FOR USE – CONVERTER ABUTMENT

The converter abutment is used in prosthetic restorations prepared by dentist in dental clinic or by technicians in dental laboratory. Converter abutments are used for changing the connection from internal to external through the Morse taper connection to secure implant-to-abutment connection. These elements are placed over the dental implant acting as part of the prosthesis and are connected to the implant by a captive screw.

After placing in implant, use adequate DESS® part in ratchet to tighten the abutment: for non-engaging abutments use octagonal key. For engaging abutments use corresponding screwdriver tip in adapter.

For both adapters torque recommended can be found in Clinical Screw Recommended Torque document; continue with standard prosthetic procedure.

APPLICATIONS AND DIRECTIONS FOR USE – TEMPORARY ABUTMENT

The temporary abutment is used for prosthetic restorations prepared by dentist at the chairside, and works as support for provisional screwed crowns, bridges and prosthesis. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend an X-ray to verify the correct settlement after tightening. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before placing the temporary restoration check the torque according to implant connection and size. It's important to seal the access hole with Teflon tape or gutta-percha and composite. Any material excess near implant platform should be removed, as it can provoke peri-implantitis and can lead to implant loss. It is very important to check static and dynamic occlusion avoiding excessive loading forces which may prolong or prevent proper osseointegration, especially in cases of immediate load.

APPLICATIONS AND DIRECTIONS FOR USE – SCAN ABUTMENT

Scan abutment is used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentist at the chairside in the CAD/CAM process. It is indicated in order to obtain geometric data from the master model using a desktop laboratory 3D scanner or for optical impressions using intraoral 3D scanner. For proper process a DESS® digital library is needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Before scanning the abutment and analogue they should be visually checked for surface damage or implant connection for any tissue residues. For greater scanning precision, we recommend to locate the flat surface of the scan abutment in palatal/lingual orientation. Fasten the abutment using corresponding screw (indications in our catalogue and on the product label) by hand or with maximum 10 Ncm. In case of intraoral scan abutment, fix it with fingers using embedded screw. Scan abutment is a precision tool and being tightened hard may change its geometry causing errors in scanning process and discrepancy in accuracy. For most scanners a scan-enhancing spray is not required.

The scanning process should be performed as recommended by CAD/CAM system manufacturer. It is important to choose the correct implant connection in software and corresponding type for the chosen restoration (engaging/non-engaging). After process the scan abutment can be loosened and placed gently on the tray or box. To complete the scanning gathering additional information (e.g. occlusion index, silicone bite, gingiva shape) can be necessary.

In case of intraoral use, sterilization is needed.

APPLICATIONS AND DIRECTIONS FOR USE – DESS LOC

DESS LOC is a resilient implant attachment used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentist at the chairside. Indicated for removable dentures.

Procedure: Measure the tissue thickness from the implant platform to the highest side of the gingival crest near the implant. Choose the corresponding height or 1.5 mm higher (the attachment should not be submerged) and using recommended torque on Clinical Screw Recommended Torque document. The retention male part can be mounted in dental clinic as well as in dental laboratory with self-curing resin using standard techniques.

APPLICATIONS AND DIRECTIONS FOR USE – PRE-MILLED ABUTMENT

Pre-milled abutment is used for prosthetic restorations prepared by dental technicians in dental laboratory. The pre-milled implant-abutment interface of the blank is fabricated on an industrial lathe in order to provide consistent and accurate results. The connection is more accurate and the abutment shape can be milled with burs to meet specific patient requirements.

For Pre-milled Abutments: The design parameters for the CAD/CAM Blank custom abutments are as follows:

Minimum wall thickness – 0.45 mm

Minimum post height – 4.0 mm


Maximum gingival height – 6.0 mm

All Blank Pre-milled abutments are for straight abutments only.

APPLICATIONS AND DIRECTIONS FOR USE – TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. By use of exchangeable screwdriver tips may be used for various implant systems. For intraoral use, all tools must be thoroughly autoclaved before use. Where using removable heads, make sure that it cannot come loose from the driver. The driver should always be fastened with a safety pin and cord attached the groove designed for this purpose. It should be handled with care, taking the necessary precautions not to hurt the patient. Before using the driver, check that it is the correct size and shape for the screw head that will be used. Only use drivers that are in perfect conditions. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

In case of adaptors please assure that corresponds to 4 x 4 connection and to its abutment system (Multi-unit®, UniAbutment®, octagonal or Locator®).

INTERNAL ANK SPECIAL REQUIREMENTS

The thickness of the anti-rotation elements is reduced due to the widened screw channel in order to accommodate the use of a normal screw. To prevent bending of those parts under applied pressure, please ensure that the screw is threaded through the abutment whilst restoration takes place and before placing it in the analog or implant. If this precaution is taken, the anti-rotation elements remain supported and safe from shear and compression forces. If, for any reason the screw needs to be removed, place it back again in its position before applying any loads on the abutment.

EXPLANATION OF THE SYMBOLS THAT APPEAR ON THE LABELS

	For single use only		Batch No.		Reference
	NON Sterile product		Read instructions	 0051	Medical product compliance marking controlled by the notified body number 0051
Rx Only			Use only on prescription		
	Manufacturer information			Manufacturer: Terrats Medical SL Calle Mogoda 75-99 08210 – Barberá del Vallés -Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 info@dess-abutments.com	

All DESS® products are manufactured under an ISO 13485-certified quality system.

3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA
 ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE® and FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants
 BioHorizons® registered trademark of BioHorizons, INC.
 Brånemark System®, NobelReplace® and Replace Select™ and NobelActive™ and NobelReplace® CC and NobelSpeedy™ and Multi-unit® registered trademarks of Nobel Biocare AB, Sweden
 ITI®, Straumann® and synOcta® registered trademarks of Straumann Holding AG, Switzerland
 MIS® registered trademark of MIS Implant Technologies Ltd
 Zimmer Screw-Vent® and SwissPlus® registered trademarks of Zimmer Dental INC, USA
 Dyna® registered trademark of Dyna Dental Engineering BV
 Locator® registered trademark of Zest Anchors INC

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. DESS® items classified as Class III devices are subject to licensing requirements as applicable in Canada. Some DESS® items may not be available in Canada - Please consult your local distributor. MetAlive-coated products are not available in Canada. These instructions replace all previous editions. For any detailed information about DESS products please contact your local distributor

Prosthetic Fixtures for dental implants

INSTRUCTIONS FOR USE

ENGLISH (UNITED STATES OF AMERICA)

WARNING – PLEASE READ CAREFULLY

U.S. Federal law restricts those devices to sale by or on the order of a licensed dentist / physician. DESS® products should only be used by dental specialists with experience in maxillary implantology and other specialties, such as dental diagnosis, planning, dental surgery or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. DESS® products designed for single-use only must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health like patient tissue infection. All DESS® components must be dry fitted before use to check the correct fitting. The clinician will be responsible for correct application of those restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please call the distributor. We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must a screw check. If the screws have suffered unusual wear, the complete integrity of the implant-abutment should be checked. New screws must be used for any revision, correction or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty. **During any intraoral use and manipulation all DESS® products must be secured to prevent aspiration due to its shape and size.**

INDICATIONS

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the laboratory.

CONTRAINDICATIONS

All materials used are biocompatible; however some patients may present allergies or hypersensitivity to any of the materials and its components (specified in the table).

COMPATIBILITY INFORMATION

All DESS® components are available for a variety of connections. For compatibility with dental implants and implant analogues please follow our catalogue or contact your local distributor.

To find necessary information about ANKYLOS® C/X system, please read the last paragraph of these instructions.

STERILIZATION AND RE-USE

All products are supplied NON STERILE. For sterilisation, we recommend autoclaving the product at 121°C for 30 minutes, drying time 30 minutes (in accordance with standard UNE-EN ISO 17665-1:2007). Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

ARTICLE/CLASS	GROUP DESS REF. Nº	MATERIAL	STERILIZATION	RE-USE
ANALOGUE / Class I	14.XXX and 14.XXX-P10	Surgical stainless steel AISI-303	In the autoclave	Not recommended
TRANSFER / Class I	17.XXX	Body: surgical stainless steel AISI-303 screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	After sterilization
SCREW / Class IIb	19.XXX and 19.XXX-P10	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
CASTABLE ABUTMENT / Class IIb (28.XXX, 29.XXX), Class I (11.XXX-P10, 12.XXX-P10)	28.XXX, 29.XXX and 11.XXX-P10, 12.XXX-P10	Polyoxymethylene copolymer (POM)	N/A	Single use
TITANIUM INTERFACE* / Class IIb	15.XXX and 16.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
CLINICAL ABUTMENTS* (STRAIGHT AND ANGLED)/ Class IIb	13.XXX, 22.XXX and 23.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
TRANSEPITHELIAL ABUTMENTS* (Multi-unit* and UniAbutment*) / Class IIb	40.XXX/X	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
HEALING ABUTMENT* / Class IIb	20.XXX/X and 20.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
CONVERTER ABUTMENT / Class IIb	34.001 and 35.001	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
TEMPORARY ABUTMENT*/Class IIb	24.XXX and 25.XXX	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use
SCAN ABUTMENT / Class I	50.XXX	Polyetheretherketone (PEEK) For ANK System: Ti alloy ELI Ti-6Al-4V	In the autoclave	In case of intraoral scanning: after sterilization
TOOLS / Class I	MDTXXX and DTXXXXXX	Handles: surgical stainless steel SUS 316L Screwdriver: stainless steel AISI 420 MOD	In the autoclave	In case of intraoral use: after sterilization

*Some parts are covered with MetAlive® coating and marked with ML letters on the beginning of the reference number.

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGUE

The analogues are for simulation of dental implant connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analogue are securely affixed. Check that the analogue and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analogue, because such a connection may be impossible to achieve with a product that may no longer be within specifications.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER

Transfer is used for closed and open tray technique, for the transfer of implant position from the intraoral situation to the model in dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove abutment, clean the connection with spray water/air and dry. Prepare transfer adequate to implant connection and platform size. Place it on implant and check correct position. Put the screw and fix manually.

- For closed tray technique choose the short screw. After tightening block the hex of the screw with wax. Use hydrocolloid, polyethylene or soft silicone. After hardening place the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For open tray technique choose the long screw. After tightening place the impression tray to ensure full access to transfer screws from outside. Before impression taking the transfers can be ferulized with acrylic resin placed over dental floss between consecutive implants. After hardening release all screws and take out the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over implant or analogue. Make sure to secure the parts with corresponding screw and observe specified torque value placed on the label or ask your local DESS® distributor. You can download torques brochure from <http://www.dess-abutments.com/downloads>. For best results, the following conditions must be meticulously met:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The driver should be placed in the longitudinal axle of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every checks thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case. If you have any doubts about the fit of the screw, consult the local distributor.
- Position the patient to avoid aspiration in case the screw falls during screwing/unscrewing.**
- Check compatibility of the screw with the implant model to which it will be connected.

Some screw models are available with DLC coating [Diamond-Like Carbon, Chromium Carbonitride coating – CrCN] to provide a low friction surface, which improves the preload of the screw and provide better hold. The Chromium Carbonitride coating is not toxic and has excellent mechanical, tribological and biological properties. DLC is often used in blood contacting implants such as stents and heart valves and to reduce wear in load bearing joints.

APPLICATIONS AND DIRECTIONS FOR USE – CASTABLE ABUTMENT

The castable abutment is for prosthetic restorations prepared by dental technicians in the dental laboratory. For greater connection precision use metal abutments which have been previously mechanised. However, if you decide to use the casting method, please take the following precautions:

- Use a flash of wax that will create sufficient gap around the abutment to compensate for the greater dilatation coefficient between the castable material and its surroundings.
- Use gentle torque to prevent deformation.
- Perform castings of a size and form that favours filling of the cavities and prevents the appearance of air bubbles.
- Use materials with high fluidity for castings in conflictive models

APPLICATIONS AND DIRECTIONS FOR USE – INTERFACE

The interface is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium interface is to support the bridge or individual zirconium dioxide elements made in CAD/CAM or manual milling machine. It can also be used as a connector between full-shape screwed crowns (from any material) and an implant, either as narrow abutment in some cases. Every titanium interface comes with a calibrated POM castable which can be burned out without residues if necessary. **Contraindications:** do not use the titanium interface for restorations with cantilever on a single implant, bruxism, insufficient space or direct metal-to-interface casting.

To use this product in digital dentistry DESS® library and scan abutment are needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Some dental CAD/CAM systems requires a double scan with diagnostical wax-up. If you need additional assistance with installation and/or proper use of the files and with use of piece itself, please contact your DESS® local distributor or find necessary information in corresponding instructions. All instructions can be downloaded from our website <http://www.dess-abutments.com/downloads>. To improve cement adhesion we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following recommendations of cement manufacturer. The ceramic surface of mesostructure in cementing zone should be sandblasted and cleaned/degreased. For secure grip the diameter of interface and its height should not be reduced e.g. by grinding.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL ABUTMENTS STRAIGHT AND ANGLED

The straight and angled abutments are used for the prosthetic restorations prepared by dental technicians in dental laboratory or by dentists at chairside in dental clinic. Straight and angled abutments form the core of prosthetic crowns and bridges or are used for their support. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts and that they are in the correct place and position for all secondary parts. We strongly recommend a radiography to verify the correct settlement after tightening. **Caution:** If customization is needed, do not grind the wall to less than 0.5 mm thick. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before cementing the final restoration, check the torque according to implant connection and size. When cementing use definitive luting cement or use temporary cement after sealing of access hole with Teflon tape or gutta-percha. Remove cement excess from edge of the crown to avoid peri-implantitis which can lead to implant loss. After cementation remember to check static and dynamic occlusion using conventional protocol.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSEPITHELIAL ABUTMENTS : Multi-unit® and UniAbutment®

Those abutments are used for prosthetic restorations prepared by dentist in dental clinic. A pre-manufactured dental implant abutment directly connected to the dental implant intended for use as an aid in prosthetic rehabilitation for multiple fixed-removable prosthesis or screw-retained. We strongly recommend to verify the final abutment seating using radiographic imaging.

DESS® offer two systems for this kind of rehabilitation: Multi-unit® system and UniAbutment®.

In case of straight Multi-unit® abutment: place appropriate abutment. Use plastic holder to facilitate the insertion. Tighten the abutment using DESS® key for Multi-unit® and manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). In case of angulated Multi-unit® abutment: place appropriate abutment using metal holder at desired position and tighten the abutment using system specific screwdriver in manual torque wrench (for recommended torque use Clinical Screw Recommended Torque document). Take an impression in a standard procedure and provisionalize the patient. If temporary prosthesis is not necessary place healing caps.

In case of UniAbutment®: place adequate abutment choosing different type (20° or 45° - depends of implant divergence) and height. Use attached holder screw to facilitate the insertion. Tighten the abutment using DESS® key for UniAbutment® and take impression using standard techniques (for recommended torque use Clinical Screw Recommended Torque document).

APPLICATIONS AND DIRECTIONS FOR USE – HEALING ABUTMENT

The healing abutment is used on the patient as an auxiliary transmucosal abutment. It is placed on the implant fixture before prosthetic restoration to facilitate the formation of a soft tissue sulcus. Choose the correct height to ensure its correct function and prevent the transmission of chewing forces. Before placing ensure that implant platform stays free from any tissue residues. Use gentle manual torque when installing.

APPLICATIONS AND DIRECTIONS FOR USE – CONVERTER ABUTMENT

The converter abutment is used in prosthetic restorations prepared by dentist in dental clinic or by technicians in dental laboratory. Converter abutments are used for changing the connection from internal to external through the Morse taper connection to secure implant-to-abutment connection. These elements are placed over the dental implant acting as part of the prosthesis and are connected to the implant by a captive screw.

After placing in implant, use adequate DESS® part in ratchet to tighten the abutment: for non-engaging abutments use octagonal key. For engaging abutments use corresponding screwdriver tip in adapter.

For both adapters torque recommended can be found in Clinical Screw Recommended Torque document, continue with standard prosthetic procedure.

APPLICATIONS AND DIRECTIONS FOR USE – TEMPORARY ABUTMENT

The temporary abutment is used for prosthetic restorations prepared by dentist at the chairside, works as support for provisional screwed crowns, bridges and prosthesis. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend a radiography to verify the correct settlement after tightening. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before placing the temporary restoration check the torque according to implant connection and size. It's important to seal the access hole with Teflon tape or gutta-percha and composite. Any material excess near implant platform should be removed, as it can provoke peri-implantitis and can lead to implant loss. It is very important to check static and dynamic occlusion avoiding excessive loading forces which may prolong or prevent the proper osseointegration, especially in cases of immediate load.

APPLICATIONS AND DIRECTIONS FOR USE – SCAN ABUTMENT

Scan abutment is used for prosthetic restorations prepared by dental technicians in dental laboratory or by dentist at the chairside in the CAD/CAM process. Indicated for obtain geometric data from the master model using a desktop laboratory 3D scanner or for optical impressions using intraoral 3D scanner. For proper process a DESS® digital library is needed. To download and install one of our libraries specific to your dental CAD system please visit our website <http://www.dess-abutments.com/downloads>

Before scanning the abutment and analogue should be visually checked for surface damage or implant connection for any tissue residues. For greater scanning precision, we recommend to locate the flat surface on spherical top in palatal/lingual orientation. Fasten the abutment using corresponding screw (indications in our catalogue and on the product label) by hand or with maximum 10 Ncm. Scan abutment is a precision tool and being tightened hard may change its geometry causing errors in scanning process and discrepancy in accuracy. For most of scanners a scan-enhancing spray is not required.

The scanning process should be performed as recommended by CAD/CAM system manufacturer. It is important to choose the correct implant connection in software and corresponding type for the chosen restoration (engaging/non-engaging). After process the scan abutment can be loosened and placed gently on the tray or box. To complete the scanning gathering additional information (e.g. occlusion index, silicone bite, gingiva shape) can be necessary.

In case of intraoral use a sterilization is needed.

APPLICATIONS AND DIRECTIONS FOR USE – METALIVE® COATING


Some DESS® parts are covered with Metalive® coating which enables protein-mediated attachment of the gingival tissue to the abutment surface. However, this does not change the general indications for use of the modified abutment.

Warning: Do not treat the gingival contact surface with any bonding agents, cements, cleaning agents or disinfectants. Do not clean the surface using abrasive technologies. Plaque removal should be done using soft scaling instruments. If sandblasted in laboratory, the gingival contact area must be protected with wax. Before sterilization or for cleaning, water vapour can be used as well as acetone and absolute alcohol in ultrasonic cleaner for 5 min. After cleaning or sterilizing do not touch surface with bare fingers and keep the surface hygienic. Please avoid disinfectants such as iodine, Hydrogen peroxide, chlorhexidine, polyvidone, etc. Avoid pointing strong laser or IR-radiation to modified surface. Do not remove Metalive® surface material from gingival contact area with dental burr or by polishing instruments such as diamond rubber. Be sure that the surface is not in contact with hard sharp edges or points. **Processing:** Metalive® abutments can be sterilized in autoclave at 121° C. After sterilization have the gingival tissue wound refreshed to enhance tissue repair and attachment, also after replacement of healing abutment with the permanent Metalive® abutment. Ensure maximal initial contact with gingival tissue. Ensure minimal movement of gingival tissue during healing. Patient should use soft toothbrushes and non-abrasive toothpaste. For cleaning the interdental spaces dental floss or interdental brushes can be used.

APPLICATIONS AND DIRECTIONS FOR USE – TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. By use of exchangeable screwdriver tips may be used for various implant systems. For intraoral use, all tools must be thoroughly autoclaved before use. Where using removable heads, make sure that it cannot come loose from the driver. The driver should always be fastened with a safety pin and cord attached the groove designed for this purpose. It should be handled with care, taking the necessary precautions not to hurt the patient. Before using the driver, check that it is the correct size and shape for the screw head that will be used. Only use drivers that are in perfect conditions. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

In case of adaptors please assure that corresponds to 4 x 4 connection and to its abutment system (Multi-unit®, UniAbutment® or octagonal).

INTERNAL ANK SPECIAL REQUIREMENTS

The thickness of the anti-rotation elements is reduced due to the widened screw channel in order to accommodate the use of a normal screw. To prevent bending of those parts under applied pressure, please ensure that the screw is threaded through the abutment whilst restoration takes place and before placing it in the analog or implant. If this precaution is taken, the anti-rotation elements remain supported and safe from shear and compression forces. If, for any reason the screw needs to be removed, place it back again in its position before applying any loads on the abutment.

EXPLANATION OF THE SYMBOLS THAT APPEAR ON THE LABELS

	For single use only		Batch No.		Reference
	NON Sterile product		Read instructions		Medical product compliance marking controlled by the notified body number 0051 (class IIb)
	Rx Only	Use only on prescription			
	Manufacturer information	Manufacturer: Terrats Medical SL Calle Mogoda 75-99 08210 – Barberá del Vallés -Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 info@dess-abutments.com			

All DESS® products are manufactured according to ISO 9001 and ISO 13485 and Class IIb bare the CE marking.

3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA
 ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE® and FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants
 BioHorizons® registered trademark of BioHorizons, INC.
 Brånemark System®, NobelReplace® and Replace Select™ and NobelActive™ and NobelReplace® CC and NobelSpeedy™ and Multi-unit® registered trademarks of Nobel Biocare AB, Sweden
 ITI®, Straumann® and synOcta® registered trademarks of Straumann Holding AG, Switzerland
 MIS® registered trademark of MIS Implant Technologies Ltd
 Zimmer Screw-Vent® and SwissPlus® registered trademarks of Zimmer Dental INC, USA
 Dyna® registered trademark of Dyna Dental Engineering BV
 Locator® registered trademark of Zest Anchors INC

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. Some DESS® items may not be available in your country. These instructions replace all previous editions. For any detailed information about DESS products please contact your local distributor.

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ITI®, Straumann® et synOcta® sont des marques déposées de Straumann Holding AG, Suisse.
MIS® est une marque déposée de MIS Implant Technologies Ltd.